stryker

Stryker CMF MEDPOR®

Customized Implant

Instructions for Use

Legend

	Is used for an instructional step that must be performed.
<ol style="list-style-type: none">1.2.3.	Is used for instructional steps that must be performed in a sequence.
<ul style="list-style-type: none">•	Is used for unordered list items.
<i>Italic</i>	Is used for references and for table or figure titles.
	Is used to supplement or clarify information.
 WARNING	Indicates a hazardous situation that, if not avoided, could result in death or serious injury.
 CAUTION	Indicates a hazardous situation that, if not avoided, could result in minor or moderate injury.
NOTICE	Indicates information that is considered important, but not hazard-related, e.g. messages relating to property damage.

Table of contents

1	Introduction	1
1.1	About this document.....	1
1.2	Other applicable documents	1
1.3	Disclaimer of liability.....	1
1.4	Definitions of terms, abbreviations and symbols	2
2	Safety information	5
3	Product information	11
3.1	Intended use/indications for use	11
3.2	Contraindications.....	11
3.3	Patient target group	12
3.4	Clinical benefit.....	12
3.5	User group	12
3.6	Adverse effects	13
3.7	For use with.....	13
3.8	Material information	13
4	Instructions for use.....	14

1 Introduction

1.1 About this document

This document is the most comprehensive source of information for the safe and effective use of the product. Read this document carefully and pay special attention to safety information. Keep this document accessible to users.

1.2 Other applicable documents

In addition to this document, the following documents are provided:

- Operative technique for Stryker CMF MEDPOR Customized Implant.
- The summary of safety and clinical performance is available in the European database on medical devices: www.ec.europa.eu/tools/eudamed.

For further information on Stryker's Electronic Instructions for Use (eIFU), go to www.ifu.stryker.com. For additional information, see product catalogs, contact your Stryker sales representative, or contact Stryker customer service.

1.3 Disclaimer of liability

In no instance will Stryker be responsible for incidental or consequential damages related to the product.

1.4 Definitions of terms, abbreviations and symbols

The following table provides definitions of terms used in this document.

Definitions of terms

Term	Definition
Healthcare professional	A healthcare professional is a surgeon, nurse, professional caregiver, or a operating medical staff.

The following table provides definitions of abbreviations used in this document.

Definitions of abbreviations

Abbreviation	Definition
CMF	Craniomaxillofacial
CT	Computed tomography
EEC	European Economic Community
eIFU	Electronic Instructions for Use
EN	European Standard
EU	European Union
EUDAMED	European database on medical devices
IFU	Instructions for Use
HDPE	Porous high-density polyethylene
MR	Magnetic resonance
MRI	Magnetic resonance imaging
PA	Polyamid
PE	Polyethylen
PTH	Persistent temporal hollowing
TMJ	Temporal mandibular joint

The following table provides the symbols used in this document, on the product, and/or on the product label.

Definitions of symbols

Symbol	Definition
	General warning sign: to signify a general warning.
	Manufacturer: indicates the medical device manufacturer.
	Date of manufacture: indicates the date when the medical device was manufactured.
	Use-by date: indicates the date after which the medical device is not to be used.
	Batch code: indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalog number: indicates the manufacturer's catalog number so that the medical device can be identified.
	Serial number: indicates the manufacturer's serial number so that the medical device can be identified.
	Sterilized using ethylene oxide: indicates a medical device that has been sterilized using ethylene oxide.
	Sterilized using irradiation: indicates a medical device that has been sterilized using irradiation.
	Sterile
	Do not resterilize: indicates a medical device that is not to be resterilized.
	Do not use if package is damaged: indicates a medical device that should not be used if the package has been damaged or opened.
	Do not reuse: indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
	Keep dry: indicates a medical device that needs to be protected from moisture.
	Consult instructions for use: indicates the need for the use to consult the instructions for use.
	Quantity: indicates the number of medical devices in the packaging.
GTIN	Global Trade Item Number
	Medical Device
	Single Sterile Barrier System

Symbol	Definition
	MR Safe: indicates a product that poses no known hazards in all MRI environments.
Rx Only	Caution: federal law (USA) restricts this device to sale by or on the order of a physician.
	Conformity assessment according to Annex I of Medical Device Regulation (EU) 2017/745 or Medical Device Directive 93/42 EEC with involvement of the notified body.

2 Safety information

This section presents general product-related safety information and recommendations for the safe and effective use of the product. For specific safety information, instructions, and recommendations, refer to each section.

The user and/or patient must report any serious product-related incident to both the manufacturer and the national competent authority where the user and/or patient is established.

General

WARNING

A lack of understanding, adequate education and experience in choice and placements of implants may lead to patient harm.

- ▶ Thoroughly understand all aspects of the surgical procedure and limitations of the implant.
- ▶ To prevent clinical failure and to maximize product life, strict adherence to the indications, contraindications, precautions and warnings is essential.
- ▶ Consider patients anatomy, up to date diagnosis and CT scan, soft tissue coverage, bone stock, existing hardware, sensitive tissues, surgical access and fixation options.
- ▶ Refer to the design proposal for a detailed description and restrictions of the implant design.

The use of a product with damaged packaging, unintentionally opened packaging, damaged sterile barrier system, or with an expired shelf life date may lead to improper use of a non-sterile or damaged product. This may lead to patient harm.

- ▶ Upon initial receipt and before first use, visually inspect the packaging for damage and confirm the integrity of the sterile barrier system for products that are received sterile.
- ▶ Dispose of this product, if the product is damaged, unintentionally opened, sterile barrier system is damaged, or the end of the product life cycle has been reached.

The product is sterile and not reprocessible. Reprocessing of a sterile product may lead to product malfunction. This may lead to patient harm.

- ▶ Do not reprocess the product.

MRI**CAUTION**

The product is MR safe. If the MR compatibility of the fixation system is not taken into consideration during MR imaging, this may result in heating and migration of the fixation products. This may lead to patient harm. For safety information, recommendations for the safe and effective use, and MRI information of the compatible fixation system plates and screws, refer to the IFUs.

Before use**WARNING**

Incorrect clinical judgement or improper planning of the procedure by the healthcare professional may lead to patient harm.

- ▶ Ensure that tissue health is satisfactory, adequate tissue coverage is available, unconfined subcranial spaces, and expected swelling space can be achieved before use.
- ▶ Ensure access to the reconstruction plan, necessary products and information during surgery.

For single stage: when a defect is created or exposed intraoperatively that is outside the boundaries of the approved virtual reconstruction plan, whether due to intentional deviation or altered patient anatomy, the implant provided may not be used to reconstruct the patient anatomy as this may lead to patient harm.

- ▶ Inspect the current reconstruction plan before the surgery.
- ▶ Ensure that an alternative reconstruction solution is available.

Unavoidable microbiological transfer or daily transient bacteremia may lead to infection. Dental manipulation, endoscopic examination and other minor surgical procedures have also been associated with transient bacteremia.

- ▶ To prevent infection at the implant site, use antibiotic prophylaxis before and after each procedure.

Intraoperative use

WARNING

An incorrect implant orientation may lead to patient harm.

- ▶ Prior to resection or modification, verify the implant orientation and position considering also less distinctive contour and curvature, CT scan, and design plan.

Inappropriate product handling or fixation may lead to sensitive tissue impingement at the patient or injury of healthcare professionals.

- ▶ Pay special attention during resection, placement and fixation of the product.
- ▶ Observe the clearance to sensitive tissues.
- ▶ Carefully handle sharp edged products.
- ▶ Deburr sharp edges generated during modification of the product.
- ▶ Ensure that the tip of the screw does not protrude past the underside of the implant or patient bone.
- ▶ Ensure that the product fits freely in the defect/void.

Exposure or penetration of the sinuses may lead to patient injury or infection.

- ▶ Do not expose the sinuses during surgical site exposure.
- ▶ Do not puncture the dura during surgical site exposure.

For single stage: resection outside of the boundaries of the approved virtual reconstruction plan may result in an inappropriate fit or implant fixation failure. This may lead to patient harm.

- ▶ Do not exceed the boundaries of the approved virtual reconstruction plan.
- ▶ Mark and resect conservatively to preserve bone stock.

Excessive intraoperative handling forces to the implant may result in implant failure. This may lead to patient harm.

- ▶ Carefully hold the implant. Hold the implant in the palm of the hand to prevent scratching or damaging the implant edges while centering the implant over the defect.
- ▶ To prevent implant damage, do not force the implant into the host bone model.
- ▶ Prevent excessive bending of the implant.
- ▶ Use low instrument speeds and hand pressure during product modification.

Implantation of host bone material may lead to patient harm.

- ▶ Do not implant the host bone or parts of the host bone.
- ▶ Do not modify the host bone.

**WARNING**

Excessive weakening of the implant during modification may lead to implant crack or break or may result in short or long term implant failure.

- ▶ Do not drill suture holes.
- ▶ Maintain a minimum of 4 mm thickness throughout regions of the implant intended to protect the brain.
- ▶ At the center joint for multi-piece implants, maintain at least a 6 mm overlap of the flange.
- ▶ At the perimeter, maintain at least 4 mm of flange contact to bone.

Any modification of the product or the bone may generate particles that may lead to puncture/irritation of the dura, or other sensitive tissues, or compromise the fit of the implant. This may lead to patient harm.

- ▶ Prevent generated particles from entering the patient.
- ▶ If needed, use suction products.
- ▶ Perform any required modification of the implant away from the surgical site using the host bone model for anatomical reference.
- ▶ After modification, remove all particles from the products, gloves, and surgical site.

Usage of high speed instruments during resection and modification of the product may lead to swirling debris or broken parts. This may lead to patient and healthcare professionals harm.

- ▶ Use protective equipment when using burring/drilling products with high speeds.
- ▶ Approach from low speed and do not exceed manufacturers limits.

Use of over-modified or mishandled implant may lead to suboptimal fixation, compromised cosmetic result or severe tissue irritation.

- ▶ Discard of over-modified or mishandled implants and use the sterile back-up implant.
- ▶ Ensure that an alternative reconstruction solution is available and inspected before the surgery.

**WARNING**

Fixation of the product into a pre-existing product or into a non-suitable area may lead to patient harm.

- ▶ Observe the area intended for implanting the product for pre-existing products and proper clearance to sensitive tissues.
- ▶ Ensure that the product fits freely inside the defect/void.
- ▶ Do not fixate the product into a pre-existing product.

A suboptimal fixation may lead to implant loosening, severe tissue irritation and necessity of a revision surgery.

- ▶ Use only compatible products along with Stryker CMF MEDPOR Customized Implants. For further details, refer to Section „3.7 For use with.“
- ▶ Use self-drilling screws of minimum of 4 mm in length.
- ▶ Attach the implant to healthy bone consistent with accepted standards for fixation.
- ▶ To prevent loosening of the implant, apply secure tightening of the fixation screws.
- ▶ To prevent tearing of the flange, center the screw and pay special attention to not over torque the screw.
- ▶ For flanged cranial implant, use at least three points of fixation, positioned in such a way that each screw is centered in the flange, maintaining at least 4 mm of material overlap to the patient bone.
- ▶ For inlay cranial implant, use at least four points of fixation distributed evenly across the implant perimeter as anatomy permits.
- ▶ For augment implant, use at least two points of fixation.

**CAUTION**

Modification of the implant with a burr may lead to collapsing or closure of pores in the area modified, potentially reducing soft tissue ingrowth over time in this section of the implant.

- ▶ Approach from low speed and do not exceed manufacturers limits.

After use**WARNING**

Patient not following the healthcare professionals' instructions may compromise the implant's fixation and result in the implant loosening or fracturing. This may lead to patient harm.

- ▶ Healthcare professionals must direct the patient to follow-up care and treatment. Healthcare professionals must instruct the patient to limit activities and to protect the implant from unreasonable strain. The patient must wear protective head gear for all sports or activities that involve possible head trauma.
- ▶ Pay special attention when bone healing has not occurred in trauma cases.
- ▶ Healthcare professionals must advise the patient to report any unusual changes of the operated site to their surgeon. The patient must be immediately evaluated by a physician if a change at the fixation site has been detected. The surgeon must evaluate the possibility of clinical failure and discuss the need for any measures deemed necessary to aid healing.

Reimplanting of the product may lead to patient harm. Forces exerted by explantation and reimplanting may compromise the physical integrity and/or dimensions of the products. Sterility cannot be assured for reused products.

- ▶ Use a new implant or alternative reconstruction solution. Make sure, that it is available and inspected before the operation.

In an unplanned explantation, inappropriate dissection of ingrown tissue may lead to patient harm.

- ▶ Pay special attention to dissection of the proximal site of the implant for example subdural site.

Disposal**WARNING**

Incorrect disposal of a product may lead to healthcare professional harm or to environmental damage. A sharp product can cause injuries to the surgical team or to the cleaning team. This may result in the spread of blood-borne pathogens.

- ▶ Dispose of products in accordance with local regulations such as the European waste catalog EAK 180103 (Cat. IB, [9, 10]).
- ▶ Pay attention to the sharp edges of products to prevent injuries and cross-contamination.
- ▶ Dispose of sharp products and empty glass items in a sharps container.
- ▶ Use specially marked lockable containers and (if applicable) protective caps.

3 Product information

The Stryker CMF MEDPOR Customized Implant is a customized patient-specific implant for cranial and/or facial application based on CT-data and input by the surgeon. Surgeons supply Stryker with a patient CT scan and the implant is designed to fit the exact patient needs. The implant is designed to fulfill surgeon's requirements and is approved by the requesting surgeon.

Each Stryker CMF MEDPOR Customized Implant Kit contains two identical sterile MEDPOR implants, and, if requested, one sterile host bone model with the defect area. The optional host bone may serve as a model in the surgical field to test the orientation and fit of the implant(s) before implantation.

The MEDPOR implants are sterilized using ethylene oxide. The host bone model is sterilized using irradiation.

The Stryker CMF MEDPOR Customized Implant is designed to maintain the interconnecting open pore structure of the MEDPOR biomaterial to allow for tissue ingrowth.

3.1 Intended use/indications for use

The Stryker CMF MEDPOR Customized Implant is indicated for the augmentation and/or restoration of bony and/or soft tissue deformities in post-traumatic, post-surgical, or congenital craniofacial defects, including, but not limited to, the correction and prevention of persistent temporal hollowing (PTH).

3.2 Contraindications

- Any active or suspected latent infection
- Any mental or neuromuscular disorder that would create an unacceptable risk of prosthesis instability, prosthesis fixation failure, or complications in postoperative care
- Bone stock compromised by poor blood supply, disease, infection or prior implantation, which cannot provide adequate support and/or fixation to the implant
- Patients with plastic/polymer allergies and/or foreign body sensitivity
- Fractures of a severely atrophic bone
- The Stryker CMF MEDPOR Customized Implant should not be used in load bearing areas such as joints like the TMJ.

3.3 Patient target group

Adult patients identified to undergo surgical procedure for the augmentation and/or restoration of bony and/or soft tissue deformities in post-traumatic, post-surgical, or congenital craniofacial defects.

3.4 Clinical benefit

The potential clinical benefits for MEDPOR Customized Implants for cranial and facial use include:

- Augmentation and/or reconstruction of post-traumatic, post-surgical or congenital defects of the cranial and/or facial skeleton and/or soft tissue with low overall complications
- Correction or prevention of persistent temporal hollowing
- Optimization of the surgical procedure and reduction of patient risks through a single-stage design
- Stable, safe and long lasting surgical results
- Good cosmetic and aesthetic satisfaction

3.5 User group

These products are intended to be used only by specialized healthcare professionals. The healthcare professionals must be licensed to perform surgery in the respective field of medicine and must be familiar with the principles of surgical procedure. The healthcare professionals must prepare the products depending on the surgery being performed and must familiarize themselves with the instructions for use in order to ensure the proper operation of this product.

The healthcare professional performing any procedure is responsible for determining the appropriateness of using the product and the specific techniques to utilize for each patient.

3.6 Adverse effects

- Patient dissatisfaction with cosmetic result
- Irritation to peri-implant sensitive tissue such as nerve, eye, dura or damage due to surgical trauma
- Immunological or foreign body response
- Early or late infection; both deep and/or superficial
- Peri-implant tissue damage during unplanned explantation
- In many instances, adverse results may be clinically related rather than implant related

3.7 For use with

The Stryker CMF MEDPOR Customized Implant products are produced and designed to be used with the standard Universal Neuro, Universal Midface, and Universal Upperface portfolio of self-drilling screws and corresponding plates. If necessary, use sharp scalpels or Stryker's burrs/drills to modify the Stryker CMF MEDPOR Customized Implant to the desired shape.

WARNING

The use of products from other manufacturers along with a Stryker product may involve incalculable risks and/or contamination of the product material and misalignments of products. This may lead to patient and healthcare professional harm.

- ▶ Do not use products from other manufacturers along with a Stryker product.

3.8 Material information

This section presents material information which applies to every product of the system unless specified otherwise.

- Stryker CMF MEDPOR Customized Implant: 100% biocompatible HDPE MEDPOR bio-material
- Host bone model: 100% polymer resin. Host bone model is optional for US only.
- Packaging material: PA/PE, cartilage, paper

4 Instructions for use

For further details on how to use the product, refer to “Operative Technique”, go to www.ifu.stryker.com.

- i** It may be necessary to remove the implant and its fixating screws from the patient in an additional surgery. If an implant removal is required, follow standard surgical procedures for explantation.

Stryker Corporation or its divisions or other affiliated entities own, use or have applied for the following trademarks or service marks: Leibinger, Stryker. All other trademarks are trademarks of their respective owners or holders.

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Copyright © 2022 Stryker
Stryker Leibinger GmbH & Co. KG