stryker

Precision Targeting System

REF 6000-801-000

REF 6000-802-000

REF 6000-803-000

REF 6000-804-032


REF 6000-805-021

REF 6000-804-060

REF 6000-806-000

Instructions for Use

Rx Only

EN		Instructions for use.....	1
DA		Brugsanvisning	28
DE		Gebrauchsanweisung.....	55
EL		Οδηγίες χρήσης.....	82
ES		Instrucciones de uso	109
FI		Käyttöohjeet	136
FR		Mode d'emploi.....	162
IT		Istruzioni per l'uso.....	189
JA		使用説明書	216
KO		사용 설명서	242
NL		Gebruiksaanwijzing	268
NO		Bruksanvisning	295
PL		Instrukcja użycia	321
PT		Instruções de utilização.....	348
RO		Instrucțiuni de utilizare	375
RU		Инструкция по применению.....	401
SV		Bruksanvisning	428
TR		Kullanım talimatı	454
ZH		使用说明	481

0








Table of Contents

1	Introduction	2
1.1	About this document.....	2
1.2	Other applicable documents	2
1.3	Definition of terms and abbreviations	3
2	Safety information	4
2.1	Safety directives.....	4
2.2	User group	6
2.3	Indications for use	6
2.4	Contraindications.....	6
3	Product overview	7
4	Preparation for use	8
4.1	Arrow alignment of the Precision Targeting Guide	8
4.2	Preparing the Tracker Sleeve.....	9
4.3	Validating the Tracker Sleeve	9
4.4	Preparing the Precision Targeting System	10
5	Intraoperative use	13
5.1	Inserting the Tracker Sleeve.....	13
5.2	Arm alignment and Guide alignment	14
5.3	Drilling.....	14
5.4	Taking a biopsy.....	16
6	Reprocessing	18
6.1	Disassembly for reprocessing	18
6.2	Reprocessing instructions.....	18
7	Disposal, transport and storage	21
7.1	Disposal	21
7.2	Transport.....	21
7.3	Storage	21
8	For use with	22
9	Technical specifications	24
10	Definition of symbols	27

1 Introduction

1.1 About this document

This document is the most comprehensive source of information for the safe and effective use of the product. Read this document carefully. Familiarization with the user documentation for the components of the guidance system before use is important. Pay special attention to safety information. Keep this document accessible to users.

The healthcare provider performing any procedure is responsible for determining the appropriateness of using the product and for the specific technique for each patient. Stryker, as a manufacturer, does not recommend a specific surgical procedure.

The following conventions are used in this document:

- The signal word **WARNING** indicates a hazardous situation that, if not avoided, could result in death or serious injury.
- The signal word **CAUTION** indicates a hazardous situation that, if not avoided, could result in minor or moderate injury.
- The signal word **NOTICE** indicates information considered important, but not hazard-related, e.g. messages relating to property damage.

1.2 Other applicable documents

In addition to this document, the following documents are provided:

- Guide for cleaning and steam-based sterilization
- User manual provided electronically with each Stryker application

1.3 Definition of terms and abbreviations

The following table provides definitions of terms and abbreviations used in this document.

Term	Definition
Accessory	Instruments used with the guidance system in order to achieve the intended use, facilitate its use or enable its functions. For example: Trackers, instruments, calibration devices, adapters, etc. Refer to the user manual for a list of accessories.
Guidance system	The computer and navigation camera system on which the surgical software applications run.
Tracker	Instrument that is used by the system to track the position of a patient, an instrument or a C-arm fluoroscope. The position is used to compute navigational information, which is then displayed on the navigation screen.

2 Safety information

2.1 Safety directives

WARNING

Before use

- Before each use, check the product for loose components, damage, malfunction, and bent or deformed parts. Do not use if these conditions exist. Failure to comply may lead to patient or medical staff injury. If service is required, contact your Stryker sales representative.
- Prior to surgery, the instrument should be checked with the guidance system to ensure they are functioning properly.

Modifications

- Do not modify this equipment without authorization of the manufacturer.

General Information

- Read and understand this information, file it in your maintenance records. Familiarization with the guidance system before its use is important. Refer to the instructions for use of the guidance system. Only trained personnel are to use this system.
- The instrument must only be used in accordance with the instructions for use contained in this manual by authorized persons who have been fully trained in their safe and effective use. The failure to follow these instructions will void your warranty.
- If the color coating on the Linear Adjustment knobs, Rotational Adjustment knobs and Locking knobs are found to be partially or fully chipped off, contact your Stryker sales representative immediately.
- The health care provider performing any procedure is responsible for determining the appropriateness of the instrument and the specific technique for each patient. Stryker, as a manufacturer, does not recommend a specific surgical procedure.

- Performing procedures with instruments other than those specified in these instructions or outside of their intended use will compromise the navigation accuracy.
- Use only Stryker approved components and accessories, unless otherwise specified. Do not modify any component or accessory. Failure to comply may result in patient and/or healthcare staff injury.
- After assembly, check that all components fit properly and securely into each other.
- Avoid subjecting the product to serious strains, such as heavy impacts. After a heavy impact, the product must be checked for defects. The product cannot be used if there are visible defects.
- ALWAYS use a sterile drape to separate a sterile device from a non-sterile zone.
- ALWAYS use a correct and rigidly installed patient tracker or a table interfaces like skull clamp.
- In case of known pre-existing conditions of variants of Creutzfeldt-Jakob-Disease (CJD) of the patient, ALWAYS ensure that the product is quarantined and not being used further. ALWAYS follow the recommendation of the national authorities (eg. WHO, RKI or CDC) concerning Creutzfeldt-Jakob-Disease.
- Do not expose the product to a high magnetic field such as from a magnetic resonance imaging (MRI) device.

 **CAUTION**

- Prior to surgery, ALWAYS ensure that all required components are available for the procedure.
- ALWAYS handle the equipment with care. DO NOT drop the device.

 The user and/or patient should report any serious product-related incident to both the manufacturer and the national competent authority of where the user and/or patient is established.

2.2 User group

Healthcare professionals (surgeon/resident, nurse/professional caregiver) educated in computer-assisted surgery and thoroughly familiar with the instructions for use and with the operation of this product.

To request an additional in-service instruction, contact Stryker.

2.3 Indications for use

United States and Canada*

Refer to the user manual supplied with the Cranial Guidance Software for system indications and contraindications.

Rest of World*

The Precision Targeting System is indicated for use as an accessory to the Cranial Guidance Software. It is intended to be used to provide guidance during drilling and biopsy of cranial tissue.

The system is indicated as an aid for locating anatomical structures in either open or percutaneous neurosurgical procedures. The system is indicated for any medical condition in which the use of computer-assisted surgery may be appropriate and where reference to a rigid anatomical structure can be identified.

*Note: Product may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Contact your Stryker representative for product availability.

2.4 Contraindications

None known.

3 Product overview

Figure 0

Part	
A	Precision Targeting Arm (6000-801-000)
A.1	Quick connect push button
A.2	Central locking knob
A.3	Starburst knob
B	Precision Targeting Guide (6000-802-000)
B.1	Precision Targeting Guide push button
B.2	Rotational adjustment knob for precise angular adjustment
B.3	Linear adjustment knob for precise linear adjustment
B.4	Locking knob for locking the trajectory
C	Tracker Sleeve (6000-803-000)
C.1	Tracker Sleeve push button
C.2	Tracker Sleeve Body (6000-803-010)
C.3	Validation Tip (6000-803-020)
D	Arm Connector Component (6000-806-000)
D.1	Starburst knob
E	Twist Drill Sleeve (6000-804-032)
E.1	Twist Drill Sleeve Body (6000-804-130)
E.2	Twist Drill Depth Stopper (6000-804-230)
F	Elite Drill Sleeve (6000-804-060)
F.1	Elite Drill Sleeve Body (6000-804-160)
F.2	Elite Drill Depth Stopper (6000-804-260)
G	Needle Sleeve (6000-805-021)

4 Preparation for use

4.1 Arrow alignment of the Precision Targeting Guide

- ✓ Prior to setting up the Precision Targeting System adjust the knobs so that the marked arrows in Figure 2 and 3 match.

1. Set all locking knobs to open.
2. Locate the adjusting knobs.


Figure 1: Location of the adjusting knobs.

3. Adjust the linear adjustment knobs so that the marked arrows in Figure 2 match. The nominal position of the arm is achieved by aligning the arrow heads.


Figure 2: Adjusting the wheel for arrow alignment.

4. Adjust the rotational adjustment knobs so that the marked arrows in Figure 3 match. The nominal position of the arm is achieved by aligning the arrow heads.


Figure 3: Adjusting the wheel for arrow alignment.

4.2 Preparing the Tracker Sleeve

1. Prepare the Tracker Sleeve by attaching the 4 Navigation Spheres to the tracker post.


Figure 4: Attaching the Navigation Spheres to the tracker post.

4.3 Validating the Tracker Sleeve

- i** For instructions on how to use the Validation Tip for the validation of the Tracker Sleeve in conjunction with the software application, refer to the user manual supplied with the software application.

1. Press the push button of the Tracker Sleeve (a).
2. Insert the Validation Tip into the Tracker Sleeve (b).

→ The Validation Tip needs to be fully inserted into the Tracker Sleeve.

Figure 5: Inserting the Validation Tip.

3. Place the Validation Tip in the Calibration Body's cone. Ensure the Validation Tip touches the bottom of the cone.


Figure 6: Validation Tip in the Calibration Body's cone.

4.4 Preparing the Precision Targeting System

WARNING

MAINTAIN SYSTEM STABILITY

- ALWAYS ensure proper assembly of the quick connect feature by listening for the audible latching of the locking mechanism.

CAUTION

- ALWAYS be careful while attaching the Precision Targeting Guide to the Precision Targeting Arm and while attaching the Tracker Sleeve to the Precision Targeting Guide to avoid a pinch point.

- i** ALWAYS ensure that the drape is not caught between the starburst connection of the skull clamp, Arm Connector Component and Precision Targeting Arm during assembly.

1. Attach the Arm Connector Component to the skull clamp.


Figure 7: Attaching the Arm Connector Component.

2. Tighten the starburst knob of the Arm Connector Component.


Figure 8: Tighten the starburst knob.

- Attach the Precision Targeting Arm directly to the skull clamp or by attaching it to the Arm Connector Component.


Figure 9: Attach the Precision Targeting Arm to the Arm Connector Component.

- Tighten the starburst knob of the Precision Targeting Arm.

i Right side of Arm Connector Component is preferred location for better maneuverability of Precision Targeting Arm.


Figure 10: Tighten the starburst knob of the Precision Targeting Arm.

- Lock the central knob of the Precision Targeting Arm.

i If the central knob has been fully loosened it may require a little force to be turned in the clockwise direction.


Figure 11 Lock the central knob of the Precision Targeting Arm.

i The other starburst of the Arm Connector Component may be used for attaching other ancillary components such as image guided surgery accessories.

6. Attach the Precision Targeting Guide to the Precision Targeting Arm by pushing the quick connect button.


Figure 12: Attach the Precision Targeting Guide to the Precision Targeting Arm.


Figure 13: Precision Targeting Guide attached.

5 Intraoperative use

WARNING

MAINTAIN SYSTEM STABILITY

- Do not press the quick connect button until required. It will disengage the Precision Targeting Guide from the Precision Targeting Arm.

i For understanding the interaction between the Precision Targeting System and the software application, refer to the user manual supplied with the application.

5.1 Inserting the Tracker Sleeve

WARNING

MAINTAIN ACCURACY

- ALWAYS ensure that the Tracker Sleeve is seated completely onto the Precision Targeting Guide before proceeding.

1. Press the push button (a).
2. Insert the Tracker Sleeve into the Precision Targeting Guide (b).


Figure 14: Inserting the Tracker Sleeve.

5.2 Arm alignment and Guide alignment

WARNING

MAINTAIN ACCURATE TRAJECTORY

- ALWAYS ensure to lock the locking knob after guide alignment.
-

CAUTION

- ALWAYS carefully handle the device when the Precision Targeting Guide is near to the skull.
-

 The Precision Targeting Guide has a color coding which is used for the Guide alignment procedure in the software application. For instructions on how to do the Guide alignment and Arm alignment, refer to the user manual supplied with the software application.

 Locking knob to be hand tighten. Do not use any tool to lock/unlock the locking knob.

5.3 Drilling

WARNING

DRILL DEPTH

- ALWAYS tighten the Drill Depth Stopper knob completely to avoid any slippage.
- ALWAYS attach the Drill Depth Stopper before inserting the drill into the sleeve.

METAL DEBRIS

- ALWAYS ensure the drill is stopped when inserting or before removing it from Drill Sleeve.

FOLLOW PLANNED TRAJECTORY

- ALWAYS ensure that the Drill Sleeve is seated completely onto Tracker Sleeve before proceeding.
 - ALWAYS ensure to use the Elite Drill Sleeve for drilling via elite attachment and the Twist Drill Sleeve for drilling via Twist drill.
-

1. Press the push button (a).
2. Insert the Drill Sleeve into the Tracker Sleeve (b).


Figure 15: Inserting the Drill Sleeve.

3. Attach the Twist Drill Depth Stopper on the Twist Drill according to the drill depth determined by the user. When using the Elite Attachment attach the Elite Drill Depth Stopper accordingly.
4. Tighten the knob of the Drill Depth Stopper.


Figure 16: Attaching the Twist Drill Depth Stopper.

5. Perform drilling.
6. Press the push button (a).
7. Remove the Drill Sleeve after drilling through the device (b).


Figure 17: Removing the Drill Sleeve.

5.4 Taking a biopsy

WARNING

FOLLOW PLANNED TRAJECTORY

- ALWAYS ensure that the Needle Sleeve is seated completely onto the Tracker Sleeve before proceeding.
- After drilling, ALWAYS ensure that the device comes back to its original trajectory before taking a sample. If not, set the trajectory again.
- ALWAYS ensure to use the Needle Sleeve for taking a biopsy sample.
- DO NOT insert the Biopsy Needle further if the user feels a drag or resistance while inserting the Biopsy Needle.

i ALWAYS ensure to tighten the knob of the Needle Sleeve before taking a sample. ALWAYS set the Needle Sleeve in unlock mode before removing the biopsy needle from its sleeve.

1. Gently loosen the knob of the Needle Sleeve (**a**) by turning it counter-clockwise. **i** Do not continue turning the knob after feeling resistance. The unlock state is achieved when the knob feels loose. Further applied force may prevent the Biopsy Needle from passing through the Needle Sleeve.
2. Press the push button (**b**).
3. Insert the Needle Sleeve into the Tracker Sleeve, which is attached to the Precision Targeting Guide (**c**).


Figure 18: Inserting the Needle Sleeve.

4. Pass the Biopsy Needle through the Needle Sleeve for biopsy sampling (**a**).
5. Lock the Needle Sleeve (**b**) to lock the outer cannula of the Biopsy Needle.
6. Take a biopsy sample.
7. Unlock the Needle Sleeve (**b**) to remove the Biopsy Needle.
8. Press the push button (**c**) to remove the Needle Sleeve from the Tracker Sleeve.


Figure 19: Inserting the Biopsy Needle.

i Refer to the instructions for use of the Biopsy Needle for instructions on how to use the Biopsy Needle.

6 Reprocessing

6.1 Disassembly for reprocessing

WARNING

- Before disassembly, ALWAYS ensure to remove all components like the Biopsy Needle from within the brain.
- During disassembly, ALWAYS hold the Precision Targeting Guide with one hand before loosening the central knob of the Precision Targeting Arm with the other hand and then move the arm away from the patient.
- Clean and sterilize the instrument before first and every use. Refer to the Guide for Cleaning and Steam-Based Sterilization.

1. Remove the sleeve by pressing the Tracker Sleeve push button.
2. Remove the Tracker Sleeve from the Precision Targeting Guide by pressing the Precision Targeting Guide push button.
3. Remove the Navigation Spheres from the Tracker Sleeve.
4. Remove the Precision Targeting Guide from the Precision Targeting Arm by pressing the quick connect button.
5. Remove the Precision Targeting Arm from the Arm Connector Component by loosening the starburst knob.
6. Remove the Arm Connector Component from the skull clamp.
7. Remove the Drill Depth Stopper from the Twist Drill/Elite attachment by unlocking the Drill Depth Stopper knob.
8. Where appropriate, put all the instruments back into the Precision Targeting System Insert Tray at their designated location.

6.2 Reprocessing instructions

This section contains only product-specific information. For general reprocessing instructions, safety notes, and reprocessing equipment, refer to the [GUIDE FOR CLEANING, DISINFECTION, AND STEAM-BASED STERILIZATION \(GCDS, TD6000005750\)](#).

 Where instructions between this document and the GCDS differ, follow this document.

The **cleaning group** of each instrument according to the GCDS is:

Part	Cleaning group
Precision Targeting Guide	4
Tracker Sleeve Body	4
Needle Sleeve	4

Part	Cleaning group
Elite Drill Sleeve	3
Twist Drill Sleeve	3
Elite Drill Depth Stopper	4
Twist Drill Depth Stopper	4
Precision Targeting Arm (refer to special reprocessing instructions in this section)	4
Validation Tip	1
Arm Connector Component	4

i These instruments can be cleaned and sterilized in the Precision Targeting System Insert Tray (REF 6000-810-000). Refer to the tray IFU for handling of the tray and placement of the instruments.

WARNING

- ALWAYS open all four locking knob and all four adjustment knobs of the Precision Targeting Guide such that the lead screws of the linear adjustment knob and threads of the rotational adjustment knobs are in maximum exposed state during reprocessing. Disassemble all components with removable parts during reprocessing.

Reprocessing instructions for the Precision Targeting Arm

Follow the GCDS with the following exceptions:

i The Precision Targeting Arm can only be cleaned in combination of manual pre-cleaning and machine washing. Manual pre-cleaning only is not sufficient.

1. Manual pre-cleaning of the device:

- Wet the wipes in the cleaning solution. Dip the wipes into the solution until saturated. After wetting the wipe, the excess liquid should be squeezed out.

i The wipes should be wet but not dripping.

- Remove heavy residues from the devices with the wipes.
- Keep the central knob of Precision Targeting Arm LOOSENED while using wet wipes and clean around articulated parts at least 3 times.
- TIGHTEN the central knob of the Precision Targeting Arm.
- Remove all residues with the cleaning tools.

- Rinse the devices with water whose quality is defined in the GCDS until the device is no longer slippery to the touch. Repeat these steps until the device is visibly clean.
2. Machine cleaning and disinfection of the device:
- Place the instrument in the Precision Targeting System Insert Tray and continue as instructed in the GCDS Section 4.3 Machine Cleaning and Disinfection with a Washer-Disinfector.

 **WARNING**

- ALWAYS ensure that the central knob of the Precision Targeting Arm is TIGHTENED during machine-cleaning, disinfection, and when exposed to running water.
 - DO NOT immerse the Precision Targeting Arm in liquids.
-

3. Sterilization

 **WARNING**

- ALWAYS ensure that the central knob of the Precision Targeting Arm is LOOSENED during sterilization
-

7 Disposal, transport and storage

7.1 Disposal

Products that have been in contact with material of human origin may be infectious. Dispose of with the necessary precautionary measures in accordance with local regulations. Ensure that infected products are decontaminated prior to recycling.

7.2 Transport

WARNING

- ALWAYS transport wrapped equipment with care to prevent damaging the sterile barrier.
-

7.3 Storage

WARNING

- ALWAYS store wrapped, processed equipment in a controlled environment that avoids extremes in temperature and moisture. Refer to the Specifications section.
 - Excessive handling of wrapped equipment will increase the likelihood of damaging the sterile barrier and may lead to contamination.
-

CAUTION

- To ensure the longevity, performance and safety of this equipment, use of the original packaging material is recommended when storing or transporting this equipment.
-

8 For use with

WARNING

Use only the Stryker-approved products specified in this section.

i Refer to the software user manual for a complete list of compatible products on system level.

Description		Catalog number
Motors		
π drive		5407-100-000
π drive+		5407-300-000
Attachments		
Elite 17 cm straight attachment		5407-120-485
Elite 20 cm straight attachment		5407-120-490
Burs		
3.0 mm acorn		5820-030-030
4.0 mm acorn		5820-030-040
5.0 mm acorn		5820-030-050
5.0 mm precision acorn		5820-030-550
3.0 mm neuro (match head) soft touch		5820-107-030
3.0 mm neuro (match head)		5820-107-430
3.0 mm neuro (match head) soft touch		5820-107-030s1
3.0 mm precision neuro (match head)		5820-107-530
3.0 mm neuro (match head) carbide		5820-107-530C
3.0 mm neuro (match head) less aggressive		5820-107-430s1
4.0 mm neuro (match head) soft touch		5820-107-040
4.0 mm precision neuro (match head)		5820-107-440

Description

Catalog number

Drill

3.2 mm twist drill

-

-

Software


Cranial Guidance Software

-

6000-670-000

Passive Markers

Navigation Spheres

8000-100-004

Guidance system

Stryker Q Guidance System

8900-100-000

Instruments

Calibration Body

8000-060-003


Navigated Biopsy Needle

6000-805-800

Tray

Precision Targeting System
Insert Tray

6000-810-000

9 Technical specifications

Specifications listed are approximate and may vary slightly from unit to unit.

Materials, dimensions, and weight

Precision Targeting Guide

Materials	Stainless steel
Length	148.850 mm
Width	149.950 mm
Height	59.68 mm
Weight	718.1 g

Precision Targeting Arm

Materials	Stainless steel, aluminum.
Length	565 mm
Width	99.7 mm
Height	80 mm
Weight	2350 g

Arm Connector Component

Materials	Stainless steel
Length	159 mm
Width	62 mm
Height	60 mm
Weight	623.4 g

Needle Sleeve

Materials	Stainless steel
Length	31.3 mm
Width	20 mm
Height	64 mm
Weight	35.8 g

Elite Drill Sleeve

Materials	Stainless steel
Diameter	17 mm
Height	59.6 mm
Weight	28.1 g

Materials, dimensions, and weight

Elite DDS

Materials	Stainless steel
Length	31.3 mm
Width	20 mm
Height	10 mm
Weight	17.7 g

Twist Drill Sleeve

Materials	Stainless steel
Diameter	17 mm
Height	59.6 mm
Weight	37.1 g

Twist DDS

Materials	Stainless steel
Length	31.3 mm
Width	20 mm
Height	10 mm
Weight	19.1 g

Tracker Sleeve




Materials	Stainless steel
Length	54 mm
Width	96.5 mm
Height	73 mm
Weight	128.3 g

Validation Tip



Materials	Stainless steel
Diameter	17 mm
Height	203 mm
Weight	67.3 g

Environmental specifications



Operation

Storage

Transportation

10 Definition of symbols

The following table defines the symbols used in this document, on the product, and on the product label.

EN ISO 7010: Graphical symbols — Safety colors and safety signs — Registered safety signs

Symbol/number	Name: Definition
 W001	General warning sign: To signify a general warning.

EN ISO 15223-1: Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements

Symbol/number	Name: Definition
 5.1.1	Manufacturer: Indicates the medical device manufacturer.
 5.1.3	Date of manufacture: Indicates the date when the medical device was manufactured.
 5.4.3	Consult instructions for use: Indicates the need for the user to consult the instructions for use.
 5.4.4	Caution: Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
 5.1.5	Batch code: Indicates the manufacturer's batch code so that the batch or lot can be identified.
 5.1.6	Catalog number: Indicates the manufacturer's catalog number so that the medical device can be identified.
 5.2.7	Non-sterile: Indicates a medical device that has not been subjected to a sterilization process.
 5.3.7	Temperature limit: Indicates the temperature limits to which the medical device can be safely exposed.
 5.3.8	Humidity limitation: Indicates the range of humidity to which the medical device can be safely exposed.
 5.3.9	Atmospheric pressure limitation: Indicates the range of atmospheric pressure to which the medical device can be safely exposed.
 5.7.7	Medical device: Indicates the item is a medical device.

21 CFR 801.109

Symbol	Name: Definition
Rx Only	Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

Product-specific symbols

Symbol	Name: Definition
GTIN	Global Trade Item Number.

ASTM F2503: Standard practice for marking medical devices and other items for safety in the MR environment

Symbol	Name: Definition
	MR Unsafe: An item that poses unacceptable risks to the patient, medical staff, or other persons within the MR environment.

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Stryker Corporation or its divisions or other affiliated entities own, use or have applied for the following trademarks or service marks: Stryker, Leibinger, Q Guidance, π drive. All other trademarks are trademarks of their respective owners or holders.
Copyright © 2024 Stryker
Stryker Leibinger GmbH & Co. KG