Registration Pointer Electromagnetic

REF 8000-050-003

Instructions For Use

Rx Only

TD8000050701 Rev AE

2024-01-09

EN Instructions for use	1
DA Brugsanvisning	32
DE Gebrauchsanweisung	63
EL Οδηγίες χρήσης	94
ES Instrucciones de uso	125
FI Käyttöohjeet	156
FR Notice d'utilisation	187
IT Istruzioni per l'uso	218
JA 使用説明書	249
KO 사용 설명서	280
NL Gebruiksaanwijzing	311
NO Bruksanvisning	342
PL Instrukcja użycia	373
PT Instruções de utilização	404
RO Instrucțiuni de utilizare	435
RU Инструкции по применению	466
SV Bruksanvisning	497
TR Kullanım talimatları	528
ZH 使用说明	559

1. How to Use this Document

1.1. About this Document

This manual is the most comprehensive source of information for the safe, effective, and compliant use and/or maintenance of the product. Read and understand this manual as well as the instructions for use supplied with the respective software application before using the product or any component compatible with the product. When combined with other medical devices, the instructions for use of these devices are to be considered as well. This manual is a permanent part of the product. Keep this manual for future reference.

The following conventions are used in this document:

The signal word **WARNING** highlights a safety-related issue. Comply with this information to prevent patient and medical staff injury.

The signal word **CAUTION** highlights a product reliability issue. Comply with this information to prevent product damage.

 Supplements or clarifies information.

The following table provides definitions of abbreviations used in this document.

Abbreviation	Definition
CDC	Centers for Disease Control and Prevention
CJD	Creutzfeldt Jakob Disease
RKI	Robert Koch Institute
WHO	World Health Organization

1.2. Symbol Definitions

EN ISO 7010: Graphical symbols – Safety colors and safety signs – Registered safety signs

Symbol

Definition

General warning sign (W001): To signify a general warning.

Refer to instruction manual/booklet (M002): To signify that the user instruction manual/booklet must be read.

ISO 15223-1: Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements

Symbol	Definition
	Manufacturer (5.1.1): Indicates the medical device manufacturer.
	Date of manufacture (5.1.3): Indicates the date when the medical device was manufactured.
	Batch code (5.1.5): Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalog number (5.1.6): Indicates the manufacturer's catalog number so that the medical device can be identified.
	Serial number (5.1.7): Indicates the manufacturer's serial number so that a specific medical device can be identified.
	Non-Sterile (5.2.7): Indicates a medical device that has not been subjected to a sterilization process.

Symbol

Definition

Keep away from sunlight (5.3.2): Indicates a medical device that needs protection from light sources.

Keep dry (5.3.4): Indicates a medical device that needs to be protected from moisture.

Temperature limit (5.3.7): Indicates the temperature limits to which the medical device can be safely exposed.

Humidity limitation (5.3.8): Indicates the range of humidity to which the medical device can be safely exposed.

Atmospheric pressure limitation (5.3.9): Indicates the range of atmospheric pressure to which the medical device can be safely exposed.

Consult instructions for use or consult electronic instructions for use (5.4.3): Indicates the need for the user to consult the instructions for use.

Symbol	Definition
--------	------------

Medical device (5.7.7): Indicates the item is a medical device.

IEC 60417 Graphical symbols for use on equipment

Symbol	Definition
--------	------------

Type BF applied part (5333): To identify a type BF applied part complying with IEC 60601-1.

Product-specific symbols

Symbol	Definition
--------	------------

Quantity: Indicates the number of products in the packaging

Note: Used to supplement or clarify information

GTIN

Global Trade Item Number

21 CFR 801.109

Symbol

Definition

Rx Only

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

Directive 2012/19/EU on waste electrical and electronic equipment (WEEE)

Symbol

Definition

Indicates that the product must be collected separately and must not be disposed of as unsorted municipal waste.

Regulatory marks and logos

Symbol

Definition

Indicates conformity to Annex I of Medical Device Regulation (EU) 2017/745 without involvement of the notified body.

2. Safety Information

WARNING

- The product is delivered in a non-sterile condition. Before first use, as well as before each following use, the product must be processed according to a validated procedure.
- Unauthorized modifications of the product are forbidden for safety reasons.
- Due to its composition, the product is not to be used together with other magnetically sensitive medical products, devices or instruments (e.g. MRI).
- Before using a medical product, assure yourself of its functional safety and proper condition via a visual inspection. Do not use the product if any defects, such as breaks, cracks, deformations or wear, are detected. Particularly important parts such as tips, notches and all moveable components should be checked with extra care.
- The product is a high-precision instrument. Avoid subjecting it to serious strains, such as heavy impacts. The product cannot be used if there are visible defects. After a heavy impact, the product must be checked for defects.
- The healthcare provider performing any procedure is responsible for determining the appropriateness of using the product and for the specific technique for each patient. Stryker, as a manufacturer, does not recommend a specific surgical procedure.
- Do not use the system during the discharge of a defibrillator. All navigated instruments must be removed from the patient and the operating area beforehand.

- i** The user or patient should report any serious product-related incident to both the manufacturer and the national competent authority where the user or patient is established.

3. Product Information

WARNING

- The product may only be used for its intended purpose and in accordance with this document and all current versions of relevant system and software application documentation.
 - This document is part of the product and must be accessible to personnel at all times. It must be provided to subsequent owners or users.
-

3.1. Stryker ENT Navigation System

3.1.1. United States / Rest of World*

Intended Use

The Registration Pointer Electromagnetic is an accessory to the Electromagnetic Navigation Unit and is intended for surface registration of the patient.

Indications for Use

The Stryker ENT Navigation System is indicated for any medical condition in which the use of stereotactic surgery may be appropriate, and where reference to a rigid anatomical structure in the field of ENT surgery, such as the paranasal sinuses, mastoid anatomy, can be identified relative to a CT- or MR-based model of the anatomy.

Example procedures include, but are not limited to the following ENT procedures:

- Transsphenoidal access procedures
- Intranasal procedures
- Sinus procedures, such as maxillary antrostomies, ethmoidectomies, sphenoidotomies/ sphenoid explorations, turbinate resections, and frontal sinusotomies
- ENT-related anterior skull base procedures

Contraindications

None known.

3.2. Cranial Guidance System

3.2.1. United States and Canada*

Refer to the user manual supplied with the Cranial Guidance Software (REF 6000-670-000) for system indications and contraindications.

Rest of World*

The Registration Pointer Electromagnetic is indicated for use as an accessory with adult and pediatric patients when using the Cranial Guidance Software using electromagnetic navigation. It is intended to be used for surface registration of patient anatomy.

The system is indicated as an aid for locating anatomical structures in either open or percutaneous neurosurgical procedures where the use of stereotactic surgery may be appropriate. The system can be used for intraoperative guidance where reference to a rigid anatomical structure, such as the skull, can be identified.

Contraindications

None known.

*Note: Product may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Contact your Stryker representative for product availability.

3.3. User Group

Healthcare professionals (surgeon/resident, nurse/professional caregiver) educated in computer-assisted surgery and thoroughly familiar with the instructions for use and with the operation of this product.

To request an additional in-service instruction, contact Stryker.

4. Product Overview


Figure 1: Registration Pointer Electromagnetic

For instructions on how to use the product for a surgical procedure, refer to the instructions for use supplied with the respective software application.

5. For Use With

 **WARNING**

- Use only Stryker-approved products, unless otherwise specified.
-

For information related to compatible software applications, refer to the instructions for use supplied with the respective software application.

6. Reprocessing

Necessary steps in the reprocessing of medical products include preparation, cleaning, disinfection, testing, maintenance, inspection, and storage. Details on how to carry out these steps for the listed products will be provided below. Also consult the CDC Guidance Document Guideline for Disinfection and Sterilization in Healthcare Facilities, 2008. Also refer to the following for the purpose of additional education but not in lieu of validated reprocessing instructions: professional organizations' clinical practice guidelines or clinical guidelines of the CDC.

-  The following cleaning and disinfection process has been validated by Stryker. If other reprocessing methods are used, make sure the reprocessing cycle is validated.

- Reprocess devices as soon as practical, typically within two hours after use.
 - Reprocess devices before you return devices for service.
-

 **WARNING**

- The product requires at minimum an intermediate level of disinfection for safe use.
 - The product is not intended for sterilization.
 - Instruments used to treat patients whose diagnosis is unclear, shall be regarded as potentially contaminated with the CJD, such instruments shall be quarantined until a non-prion diagnosis is identified.
 - If the product is used on a patient who suffers from a disease whose pathogens cannot be eliminated with conventional, standardized procedures, the product must be disposed of or reprocessed following the requirements of the responsible public authority (e.g. WHO, RKI, CDC).
 - When preparing and applying chemical solutions, the manufacturer's instructions regarding mixing ratios and dwell times must be followed closely. Using false mixing ratios or extended dwell times can lead to product contamination.
 - Do not overload the washer-disinfector. Avoid unwashed areas. Ensure that the product is placed securely in the machine.
-

 **CAUTION**

- In regard to the composition of materials, ensure that the reprocessing is monomaterial. Instruments made out of titanium may change their surface color. Base metals may develop surface rust.
-

Personal Safety

Use the following personal protective equipment:

- (single-use) shirt
- (single-use) gloves
- safety glasses
- anti-splash mask

Refer to the user instructions of the cleaning agent and disinfectant for additional protective equipment.

6.1. Preparation for Reprocessing

To prevent residues from drying out, the following steps are recommended as soon as practical after use:

1. Wipe off the instruments.
-

WARNING

- Do not use any fixing agents or hot water, as this could interfere with the success of the reprocessing procedure.
-

CAUTION

- Do not use any NaCl solution for pre-cleaning, as chloride ions can lead to corrosion on metals.
-
2. Rinse the instruments thoroughly with water (water temperature max. 45 °C).
 3. During transport to the reprocessing site, store the instruments in a closed container, to prevent residues from drying out. Cover devices with a damp cloth or enclose devices in a heat-sealed bag.

4. Check the instruments for damage and wear both before and after each application. Pay special attention to breaks, cracks, deformations, distortions, and corruptions. Check particularly important parts such as tips, notches, and all moveable components with extra care.
5. Dispose of any damaged instruments.

6.2. Cleaning and Disinfection

The cleaning and disinfection can be performed manually or automated, but automated decontamination is preferred.

6.2.1. Manual Cleaning and Disinfection

Required equipment:

- Low-lint cloth, swab, cotton swabs
- Soft cleaning brushes without metal, for surgical instruments (e.g. Medisafe MED100.25, bristle length 7 mm)
- Cleaning pistol (at least 2 bar static water pressure, e.g. Stoeckert SELECTA, REF 1769534)
- Single-use syringe (20 ml)
- Alkaline cleaner (neodisher MediClean forte, Dr. Weigert, Hamburg, Germany)

- Freshly prepared purified water, highly purified water or sterile water with less than 10_cfu/ml and 0.25_EU/ml. In the U.S., refer to AAMI TIR 34 for guidance on water quality
 - High-level disinfectant based on aldehyde (Cidex OPA, Johnson & Johnson Advanced Sterilization Products, Irvine, CA, USA)
 - Medical grade compressed air and/or 70% ethyl or isopropyl alcohol
-

 **WARNING**

If the product is intended for use in known immunocompromised patients or potentially immunocompromised patients, based on institutional procedures (e.g., high risk population served), the instruments shall be rinsed with purified water, using sterile technique when rinsing and handling. For all other devices, a sterile water rinse is recommended when practical. Otherwise, potable tap water rinse is acceptable.

 **CAUTION**

- Do not use any sharp objects such as metal brushes or other metal objects for manual cleaning.
 - Inadequate rinsing or prolonged immersion in a cleaning or disinfection bath may lead to corrosion. Consult the instruction leaflet of the respective cleaning and disinfection products regarding dwell times.
 - Do not use ultrasonic baths for cleaning the electromagnetic instruments.
-

Manual Cleaning

For manual cleaning, Stryker validated the alkaline cleaner neodisher MediClean forte (Dr. Weigert, Hamburg, Germany) with 0.5 – 3.0% in purified water.

Prepare the equipment and the devices:

1. Dilute the cleaning agent at 22–40 °C. Stir the cleaning solution to homogenize it.
2. Place the dismantled instruments in a cleaning solution (22-40 °C) for 5-10 minutes. Refer to the exposure times indicated by the detergent manufacturer.

WARNING

- Completely wet all instruments with the cleaning solution, keeping them free from air bubbles. Ensure that all cavities are accessible. If necessary, enable bubbles to escape by moving and tilting the instruments.
 - Thorough cleansing is necessary to ensure successful disinfection.
 - It is essential to pay attention to the instructions for use from the detergent and disinfectant manufacturers.
-

3. Remove remaining soil with the help of the cleaning tools described above:
 - Brush off any coarse soiling from the product surface.
 - Wipe off the cable with a soft cloth soaked in purified water.
4. Flush the instruments with a cleaning pistol (purified water, 22-35 °C) for at least 20 seconds.
5. Rinse all parts at least twice for 10 seconds each under running water (purified water, 22-35 °C), in order to remove all residues.
6. Dry the instrument carefully using a low-lint cloth or swab. Hard-to-access areas can be dried using medical grade compressed air.
7. Verify that no visible staining or cleaning solution, especially in hard-to-reach areas remains on the instrument. Remove visible staining or cleaning solution by repeating the cleaning process from the beginning

 **WARNING**

- Make sure that all product parts are dried. Take special care to ensure that areas inside the product are dried.
-

Manual Disinfection

For manual disinfection, Stryker recommends high-level chemical disinfectant Cidex OPA (Advanced Sterilization Products, REF 20391, active concentration of 0.55% orthophthalaldehyde). Refer to the Cidex OPA labeling and instructions for use for additional instructions.

i Non-sterilizable, but disinfectable instruments have a red plug.

1. Immerse the instrument completely in Cidex OPA solution for a minimum of 12 minutes at 22 °C (72 °F) to destroy all pathogenic microorganisms. Make sure that you eliminate air pockets (e.g. in the plug).
2. Remove the instrument from the Cidex OPA solution, and thoroughly rinse the instrument:
 - a. Immerse the instrument completely in a large volume (e.g. 2 gallons) of sterile-filtered water, and keep the device totally immersed for a minimum of 1 minute.
 - b. Remove the device and discard the rinse water. Do not reuse the water for rinsing or any other purpose.
3. Repeat steps 2a-b twice (for a total of 3 rinses).
4. Visually inspect the device. In case of any residues of disinfectant solution, repeat step 2 a–b.

 **WARNING**

- Three (3) separate, large volume water immersion rinses are required. Disinfectant residues may cause serious side effects.
 - When using potable water for rinsing, the user should be aware of the increased risk of recontaminating the device or medical equipment with microorganisms which may be present in potable water supplies. Water treatment systems, such as softeners or deionizers, may add microorganisms to the treated water to the extent that microbial content of the water at the point of use could exceed that of the pretreated drinking water.
 - Do not use any disinfectants that contain phenols, chlorine components, or peracetic acid without corrosion protection.
-
5. Dry the product carefully using a low-lint cloth or swab. Hard-to-access areas can be dried using medical grade compressed air. For better drying, a final rinse using a 70% isopropyl alcohol solution can be used to speed the drying process.

 **WARNING**

- A device that is not completely dried provides an ideal situation for rapid colonization of bacteria. Rapid drying will avoid possible colonization but may not result in a device free from these bacteria. A final rinse using a 70% isopropyl alcohol solution can be used to speed the drying process and reduce the numbers of any organism present as a result of rinsing with purified water.
-

6.2.2. Automated Cleaning and Thermal Disinfection

A combined cleaning and disinfection process should take place in a fully automatic washer-disinfector.

Stryker recommends the cleaning agent neodisher MediClean forte (Dr. Weigert, Hamburg, Germany) and thermal disinfection methods.

Required equipment:

- Automatic washer-disinfector compliant with ISO 15883 (series) with a validated procedure (e.g. Miele G7836 CD)
- Small part sieves for safe fixation of the device
- Alkaline cleaner
(neodisher MediClean forte, Dr. Weigert, Hamburg, Germany)

- Cleaning pistol
(at least 2 bar static water pressure, e.g. Stoeckert SELECTA, REF 1769534)
 - Purified water
-

 **WARNING**

- Only use washer-disinfectors that are suitable for the cleaning and disinfection of surgical instruments.
 - Refer to the maximum load of your washer-disinfector.
 - It is essential to pay attention to the instructions for use from the detergent manufacturers.
 - To ensure the effectiveness of the cleaning process, follow the described pre-cleaning steps prior to cleaning.
-

1. Manually pre-clean the instruments:

- a. For the removal of coarse dirt, place the instruments in cold purified water for 5-10 minutes and wipe off visible staining with a soft cloth soaked in cleaning solution.
- b. Flush the instruments with a cleaning pistol (purified water, 22-35 °C) for at least 20 seconds.
- c. Visually inspect the instruments. If the instruments are not visually clean, repeat step 1a-b.

2. Place the instruments in the washer-disinfector.
 - Safely fix all instruments using a small part sieve tray from the manufacturer of your washer-disinfector.
 - Place devices such that they do not collide during cleaning.
 - Ensure that no device is obstructed by a large device (e.g. plates cup).
 - For devices with recessed features (holes, cannulae, crevices), connect the hoses of the washer-disinfector to the recessed features.
 - Place articulated devices in the fully open position and cannulated devices in the horizontal position.
 - Arrange the instruments in such a way that no areas are left unwashed and inner and outer surfaces are reached by the rinsing medium.
 - Do not overload the washer-disinfector.
3. Set the parameters on the washer-disinfector. Note the instructions from the manufacturer of your washer-disinfector for surgical instruments.

The following cleaning cycle was validated in a Miele G7836 CD (two level rack, injection rate 5 ml/l). The used automated cleaning program had the following parameters:

- At least 2 minutes pre-cleaning with purified water (16 °C)

- At least 10 minutes cleaning at 50 °C with 0.5% solution of neodisher MediClean forte (Dr. Weigert, Hamburg, Germany) in purified water
 - At least 2 minutes rinsing with purified water (22-35 °C)
 - Thermal disinfection with purified water at 90 °C for 1 or 5 minutes or in conformance with the national requirements according the A0 value (without an additional agent).
 - Drying at 110 °C for at least 15 minutes if the machine cycle includes a drying step.
4. Unload the washer-disinfector.
 5. If the machine cycle does not include a drying step, dry the devices in an oven at 110 °C for at least 15 minutes.
 6. Inspect cannulated, flexible and articulated devices for complete dryness. If required, use medical grade compressed air to further dry.

 **WARNING**

- Make sure that all product parts are carefully dried. Take special care to ensure that areas inside the product are dried.
-

 **CAUTION**

- Do not use any rinsing agents with drying aids for the last wash cycle. These could remain on the surface causing interactions that damage the instrument.
-

6.3. Visual Inspection and Functionality Testing

1. Visually inspect the products. Pay special attention to rough surfaces and hard-to-access areas. It may be necessary to use a magnifying glass.

Ensure that the products are:

- visually clean
 - dried
 - free from all contamination, such as organic material
 - free from cleaning agent residues
2. If the product parts are not visually clean, repeat the cleaning and disinfecting cycle until this is achieved. Otherwise effective disinfection cannot be guaranteed.
 3. Check the products for proper function and against damage and wear before and after each application. Check particularly important parts such as tips, notches, and all moveable components with extra care.

Pay special attention to:

- breaks
- cracks
- deformations
- corrosion

4. Dispose of any damaged instruments.

7. Maintenance


Maintenance and repair may only be conducted by the manufacturer or authorized partners. The product components may only be sent back to the manufacturer in a cleaned and disinfected condition. Sharp or pointed components need to be sent back in a protected state.

8. Disposal

Risk of injury or infection can be avoided by safe disposal of the product components. Sharp and pointed components need to be collected and locked in a tight and break-proof container. They must be stored in such a way that they are protected from unauthorized use. Contaminated products are to be supplied to a hazardous waste site and handled in a way that contamination of third parties is excluded.

- In accordance with European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) electronic products should be collected separately for recycling. Do not dispose of as unsorted municipal waste. Contact local distributor for disposal information. Ensure infected equipment is decontaminated prior to recycling.
- The plugs used in the electronic products contain the following substance: Lead, CAS No. 7439-92-1 (1907/2006 REACH). Handling Instructions: No specific precautions are required for handling items manufactured from alloys containing lead in the supplied condition.
- The silicone glue used in these products contains the following substances: Decamethylcyclopentasiloxane, CAS No. 541-02-6 and Dodecamethylcyclohexasiloxane, CAS No. 540-97-6 (1907/2006 REACH).

9. Technical Specifications

Environmental Limitations	Operation	Storage and Transportation
Temperature	 A diagram showing a thermometer with a scale. The minimum temperature is labeled as 10°C and the maximum as 30°C.	 A diagram showing a thermometer with a scale. The minimum temperature is labeled as -10°C and the maximum as 50°C.
Relative humidity	 A circular icon containing a percentage sign and wavy lines representing water. The upper limit is labeled as 75% and the lower limit as 30%.	 An icon of an umbrella with rain falling from it, indicating protection from moisture.
Atmospheric air pressure	 A circular icon with two arrows pointing towards a central point. The upper limit is labeled as 106 kPa and the lower limit as 80 kPa.	-

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Einschränkungen in Bezug auf die Umgebungsbedingungen	Betrieb	Lagerung und Transport
Temperatur	 A diagram showing a temperature range from 10°C to 30°C. A thermometer icon is positioned between the two values, with lines connecting the values to the thermometer's scale.	 A diagram showing a temperature range from -10°C to 50°C. A thermometer icon is positioned between the two values, with lines connecting the values to the thermometer's scale.
Relative Luftfeuchte	 A diagram showing a relative humidity range from 30% to 75%. A circular icon containing a percentage sign and wavy lines is positioned between the two values, with lines connecting the values to the icon.	 An icon of an umbrella with rain falling on it, representing protection from moisture or weather.
Luftdruck	 A diagram showing an air pressure range from 80 kPa to 106 kPa. A circular icon containing a pressure gauge symbol is positioned between the two values, with lines connecting the values to the icon.	-

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
9. Technische specificaties

Omgevingslimieten	Bediening	Opslag en vervoer
Temperatuur	 A thermometer icon with a liquid level. A line points from the top of the liquid to the label "30°C". Another line points from the bottom of the liquid to the label "10°C".	 A thermometer icon with a liquid level. A line points from the top of the liquid to the label "50°C". Another line points from the bottom of the liquid to the label "-10°C".
Relatieve luchtvochtigheid	 A circular icon containing a percentage sign and wavy lines at the bottom. A line points from the top right to the label "75%". Another line points from the bottom left to the label "30%".	 A simple icon of an open umbrella with rain falling from it.
Atmosferische luchtdruk	 A circular icon containing a pressure gauge symbol with two arrows pointing towards a central dot. A line points from the top right to the label "106 kPa". Another line points from the bottom left to the label "80 kPa".	-

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This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

U.S. Patents: www.stryker.com/patents

Stryker Corporation or its divisions or other affiliated entities own, use or have applied for the following trademarks or service marks: Leibinger, Stryker. All other trademarks are trademarks of their respective owners or holders.

Copyright © 2024 Stryker

Stryker Leibinger GmbH & Co. KG
Bötzinger Straße 41
79111 Freiburg (Germany)
t: +49 761 4512 0 (Germany)
t: +1 269 323 7700 (USA)